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8. Referanslar

- 1) Na Zhu, Dingyu Zhang, Wenling Wang, Xingwang Li, Bo Yang, Jingdong Song, Xiang Zhao, Baoying Huang, Weifeng Shi, Roujian Lu, Peihua Niu, Faxian Zhan, Xuejun Ma, Dayan Wang, Wenbo Xu,Guizhen Wu, George F. Gao, Wenjie Tan,et al. A novel coronavirus from patientswith pneumonia in China, 2019, N Engl J Med. 2020 Jan 24.
- 2) Severe acute respiratory syndrome-related coronavirus: The species and its viruses – a statement of the Coronavirus Study Group. Alexander E. Gorbalenya, Susan C. Baker, Ralph S. Baric, Raoul J. de Groot, Christian Drosten, Anastasia A. Gulyaeva, Bart L. Haagmans, Chris Lauber, Andrey M Leontovich, Benjamin W. Neuman, Dmitry Penzar, Stanley Perlman, Leo L.M. Poon, Dmitry Samborskiy, Igor A. Sidorov, Isabel Sola, John Ziebuhr. Severe acute respiratory syndrome-related coronavirus: The species and its viruses – a statement of the Coronavirus Study Group. bioRxiv2020.02.07.937862; doi: <https://doi.org/10.1101/2020.02.07.937862>
- 3) Laboratory testing for Middle East Respiratory Syndrome coronavirus, interim guidance (revised), January 2019,WHO/MERS/LAB/15.1/Rev1/2019, World Health Organization, 2018. (<https://apps.who.int/iris/bitstream/handle/10665/259952/WHO-MERS-LAB-15.1-Rev1-2018-eng.pdf;jsessionid=08939A780A5A4552E C8279D0E6D650E4?sequence=1>)
- 4) Managing epidemics, key facts about major deadly diseases. Geneva: World Health Organization; 2018. (<https://apps.who.int/iris/handle/10665/272442>).
- 5) WHO Global Influenza Surveillance Network Manual for the laboratory diagnosis and virological surveillance of influenza, WHO, 2011 (https://www.who.int/influenza/gisrs_laboratory/manual_diagnosis_surveillance_influenza/en/).
- 6) Protocol to investigate non-seasonal influenza and other emerging acute respiratory diseases. Geneva: World Health Organization;2018. (https://www.who.int/influenza/resources/publications/outbreak_investigation_protocol/en/).
- 7) WHO Recommended Surveillance Standards WHO/CDS/CSR/ISR/99.2 (<https://www.who.int/csr/resources/publications/surveillance/whocdscsr992.pdf>).
- 8) Guideline for the collection of clinical specimens during field investigation of outbreaks WHO/CDS/CSR/EDC/200.4 (https://www.who.int/ihr/publications/WHO_CDS_CSR_EDC_200_4/en/).
- 9) WHO laboratory biosafety manual, third edition. Geneva: World Health Organization; 2004. (http://www.who.int/csr/resources/publications/biosafety/WHO_CDS_CSR_LYO_2004_11/en/).
- 10) Global Surveillance for human infection with coronavirus disease (COVID-2019), Interim guidance, Geneva, World Health Organization, 2020. ([https://www.who.int/publicationsdetail/global-surveillance-for-humaninfection-with-novel-coronavirus-\(2019-ncov\)](https://www.who.int/publicationsdetail/global-surveillance-for-humaninfection-with-novel-coronavirus-(2019-ncov))).
- 11) Qun Li, Xuhua Guan, Peng Wu, et al. Early Transmission Dynamics in Wuhan, China,of Novel Coronavirus–Infected Pneumonia.N Eng Jour Med.January 29, 2020.

- 12) Wei-jie Guan, Zheng-yi Ni, Yu Hu, Wen-hua Liang, Chun-quan Ou, Jian-xing He, Lei Liu, Hong Shan, Chun-liang Lei, David S.C. Hui, Bin Du, Lan-juan Li, Guang Zeng, Kwok-Yung Yuen, Ru-chong Chen, Chunli Tang, Tao Wang, Ping-yan Chen, Jie Xiang, Shi-yue Li 1, Jin-lin Wang, Zi-jing Liang, Yi-xiang Peng, Li Wei, Yong Liu, Ya-hua Hu, Peng Peng, Jian-ming Wang, Ji-yang Liu, M.D, Zhong Chen, Gang Li, Zhi-jian Zheng, Shao-qin Qiu, Jie Luo, Chang-jiang Ye, Shao-yong Zhu, Nan-shan Zhong, Clinical characteristics of 2019 novel coronavirus infection in China. medRxiv preprint doi: <https://doi.org/10.1101/2020.02.06.20020974>.
- 13) Nathalie M Linton, Tetsuro Kobayashi, Yichi Yang, Katsuma Hayashi, Andrei R. Akhmetzhanov, Sung-mok Jung, Baoyin Yuan, Ryo Kinoshita, Hiroshi Nishiura. Incubation Period and Other Epidemiological Characteristics of 2019 Novel Coronavirus Infections with Right Truncation: A Statistical Analysis of Publicly Available Case Data. *J. Clin. Med.* 2020, 9(2), 538; <https://doi.org/10.3390/jcm9020538>.
- 14) XU Kaijin, CAI Hongliu, SHEN Yihong, et al. Management of Corona Virus Disease-19 (COVID-19): the Zhejiang Experience. *Jour Zhejiang Univ.* DOI : 10.3785/j.issn.1008-9292.2020.02.02.
- 15) Wei Zhang, Rong-Hui Du, Bei Li, Xiao-Shuang Zheng, Xing-Lou Yang, Ben Hu, Yan-Yi Wang, Geng-Fu Xiao, Bing Yan, Zheng-Li Shi & Peng Zhou (2020) Molecular and serological investigation of 2019-nCoV infected patients: implication of multiple shedding routes, *Emerging Microbes & Infections*, 9:1, 386-389 (<https://www.tandfonline.com/doi/full/10.1080/22221751.2020.1729071>).
- 16) Yong Zhang, Cao Chen, Shuangli Zhu, Chang Shu, Dongyan Wang, Jingdong Song, Yang Song, Wei Zhen, Zijian Feng, Guizhen Wu, Jun Xu, Wenbo Xu. Isolation of 2019-nCoV from a Stool Specimen of a Laboratory-Confirmed Case of the Coronavirus Disease 2019 (COVID-19), 2020, *CCDC Weekly / Vol. 2/ No. 8*.
- 17) Infection prevention and control during health care when novel coronavirus (nCoV) infection is suspected, interim guidance, January 2020. Geneva: World Health Organization; 2020. ([https://www.who.int/publicationsdetail/infection-prevention-and-control-during-health-care-when-novel-coronavirus-\(ncov\)-infection-is-suspected-20200125](https://www.who.int/publicationsdetail/infection-prevention-and-control-during-health-care-when-novel-coronavirus-(ncov)-infection-is-suspected-20200125)).
- 18) Laboratory biosafety guidance related to the novel coronavirus (2019-nCoV), World Health Organization; 2020 (https://www.who.int/docs/defaultsource/coronaviruse/laboratory-biosafetynovel-coronavirus-version-1-1.pdf?sfvrsn=912a9847_2).
- 19) Shi X, Gong E, Gao D, et al. Severe acute respiratory syndrome associated coronavirus is detected in intestinal tissues of fatal cases. *Am J Gastroenterol.* 2005;100 (1):169–176.
- 20) Zhou J, Li C, Zhao G, et al. Human intestinal tract serves as an alternative infection route for Middle East respiratory syndrome coronavirus. *Sci Adv.* 2017;3(11).
- 21) Ding Y, He L, Zhang Q, et al. Organ distribution of severe acute respiratory syndrome (SARS) associated coronavirus (SARS-CoV) in SARS patients: implications for pathogenesis and virus transmission pathways. *J. Pathol.* 2004;203(2):622–630.
- 22) Guidance on regulations for the transport of infectious substances 2019–2020. Geneva: World Health Organization; 2019. (<https://www.who.int/ihr/publications/WHO-WHE-CPI-2019.20/en/>).
- 23) Guidance to minimize risks for facilities collecting, handling or storing materials potentially infectious for polioviruses (PIM Guidance). Geneva: World Health Organization; 2018. <http://polioeradication.org/wpcontent/uploads/2016/07/PIM-guidance-20190122-EN.pdf>.
- 24) Meyer B, Drosten C, Müller MA. Serological assays for emerging coronaviruses: challenges and pitfalls. *Virus Res.* 2014 Dec 19;194:175-83.
- 25) Bai Shaoli, Wang Jianyun, Zhou (Yingquan, Yu Desheng, Gao Xiaomin, Li Lingling, Yang Fan. Analysis of the first family epidemic situation of new coronavirus pneumonia in Gansu Province. *Chinese Journal of Preventive medicine*, 2020, 54.
- 26) Shu-Yuan Xiao, Yingjie Wu, Juan Li, Evolving status of the 2019 novel coronavirus infections: proposal of conventional serologic assays for disease diagnostics and infection.
- 27) GISAIID.org (<https://www.gisaid.org/>), accessed on 19 February 2020

- 28) World Health Organization. International Health Regulations (2005), third edition. Geneva: World Health Organization; 2016
(<http://www.who.int/ihr/publications/9789241580496/en/>).
- 29) Protocol to investigate non-seasonal influenza and other emerging acute respiratory diseases. Geneva: World Health Organization; 2018.
(https://www.who.int/influenza/resources/publications/outbreak_investigation_protocol/en/).
- 30) Bruce et al. JCM. 2011. Evaluation of Swabs, Transport Media, and Specimen Transport Conditions for Optimal Detection of Viruses by PCR.

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DSÖ referans numarası: [WHO/COVID-19/laboratory/2020.4](https://www.who.int/laboratory/2020.4)

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² [https://www.who.int/publications-detail/global-surveillance-for-human-infection-with-novel-coronavirus-\(2019-ncov\)](https://www.who.int/publications-detail/global-surveillance-for-human-infection-with-novel-coronavirus-(2019-ncov))

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© Dünya Sağlık Örgütü 2020. Tüm hakları saklıdır.